


DYNAMIC RETRACTION

For Adjustable Retraction / Suspension

Versa Lifter Needle



Versa Lifter Needle is fully repositionable and adjustable intra-operatively: The anchoring clip can be placed and replaced to exactly adjust the retraction angle to the needs of the surgeon.

Targeted specialties: GYN, Urology, Colorectal, General Surgery


Highlighted procedures: Ovarian Cyst, Sacro-colpopexy, Sleeve Gastrectomy, Cholecystectomy

Used for suspension of: Uterus, Ovaries, Sigmoid Colon, Peritoneum, Bowel, Stomach, Gallbladder, Bladder, Appendix, Hernia Mesh...

STATIC RETRACTION

For Fixed Retraction / Suspension

TPEA Lifter



TPEA Lifter provides a simple and effective solution for lifting or retraction during laparoscopic surgery for maximum operative field exposure.

Product Type	Product Code	Specifications	Dynamic Retraction	Static Retraction	Bariatric
Versa Lifter Needle	VN10A05	5 cm Needle, 18 cm Suture, 15 x 5 mm Elongated Silicone Device, 5 per box	✓	✓	
	VN10X05	6.5 cm Needle, 18 cm Suture, 15 x 5 mm Elongated Silicone Device, 5 per box	✓	✓	✓
TPEA Lifter	TPL5A180-B10	5 cm Needle, 18 cm Suture, 10 per box		✓	
	TPL5X180-B10	6.5 cm Needle, 18 cm Suture, 10 per box		✓	✓

Features & Benefits:

- Easy to use, rapid set-up
- Hands-free solution for optimal surgical field exposure
- Allows retraction in any desired position versus fixed port position
- Cost-effective, eliminates a trocar port and an assistant
- Fully compatible with robotic procedures: prevents instrument and robotic arm interference
- Two needle lengths to accommodate all wall thicknesses
- Reduced patient discomfort and scarring, better cosmetic result
- Shorter procedure / OR time